


VISIONARY

# MESA GRANDE

*ortho-k*



*"Keep 'em outta glasses"*


## Introducing Efficiency to Ortho-K

---

ef•fi•cient  
*adjective*

a design achieving maximum productivity with minimum wasted effort or expense

---



### Available Large Range of Parameters:

- Four or five curve options allow for design flexibility and easier fit
- Large base curve to steep curve ratio guarantees the most efficient rate of correction
- Alignment curve and steep curve are easily designed from K's and refractions
- Liberal warranty options are available to ease financial stress

MESA GRANDE  
ortho-k

10/27/2019

Account Number:	
Point of Contact Name:	
Phone Number:	
Email:	
Patient Name:	

ENTER PATIENT DATA	
<b>OD LENS</b>	<b>OS LENS</b>
Flat K (D)	Flat K (D)
Rx (sphere)	Rx (sphere)
HVID (mm)	HVID (mm)
NOTES:	

RESET

MG LENS DESIGN	
<b>OD LENS</b>	<b>OS LENS</b>
POWER	POWER
BASE CURVE	BASE CURVE
DIAMETER	DIAMETER
OZ	OZ
STEEP CURVE RADIUS	STEEP CURVE RADIUS
STEEP CURVE WIDTH	STEEP CURVE WIDTH
ALIGNMENT CURVE 1 RADIUS	ALIGNMENT CURVE 1 RADIUS
ALIGNMENT CURVE 1 WIDTH	ALIGNMENT CURVE 1 WIDTH
ALIGNMENT CURVE 2 RADIUS	ALIGNMENT CURVE 2 RADIUS
ALIGNMENT CURVE 2 WIDTH	ALIGNMENT CURVE 2 WIDTH
PERIPHERAL CURVE RADIUS	PERIPHERAL CURVE RADIUS
PERIPHERAL CURVE WIDTH	PERIPHERAL CURVE WIDTH


PRINT ORDER FORM

COPY TO CLIPBOARD

Visionary Contact Lens  
info@visionarylens.com  
800-488-2020

### Simple Fitting Process:

- Simple and fast fitting using Visionary MG design tool
- Lens design based on HVID, Kerataconus Readings and Rx



### Implementation of Most Advanced Technologies:

- Made from Boston Equalens II with a high Dk of 85 for high oxygen permeability and healthy comfortable overnight contact lens wear
- Repeatable and accurate design producibility using advanced technology DAC lathes
- Available with Tangible Hydra-PEG treatment
- Quality inspection tools (Rotlex) and controls exceed FDA requirements

